



Flexion Series

Mobile C-arm X-ray System



WDM

Beijing Wandong Medical Technology Co., Ltd.

Address: Building 3 No.9 Jiuxianqiaodong Road, Chaoyang

District, 100015, Beijing, P.R.China

Tel: +86-10-84575792/3/5/6

Fax: +86-10-84575794

Web: www.wandong.com.cn

E-mail: International@wandong.com.cn

Note: Design and specification subject to change without notice.

WDM Esaote (Suzhou) Medical Technology Co., Ltd.

Address: Building 7, Yangshan Science and Technology Industrial Park,

No.8 Jinyan Road, High-tech Zone, Suzhou, 215000, Jiangsu


P.R. China

WDM

Integrated Design

All-in-one system for easy movement in the operating room.

Smaller footprint, ideal for complex OR environments and space saving.




E-Control
Flexible Operation

SMART
High Quality Imaging

I-Dose
Low Dose Control




E-Control
Flexible Operation




Longitudinal Travel: 200 mm



Vertical Travel: 400 mm




SID:1080mm Arc Depth: 708 mm Opening: 884 mm



Orbital Movement: 150°



Rotation: $\pm 225^\circ$



27 inch Wide Display



Built-in UPS for Unrestricted Transport
Single-hand Operation



100N Release/Locking Force



100N Moving Force



Color-Coded Design

SMART High-Quality Imaging




High Power · Small Spot Size ·
Continuous Operation for Long Time

0.3/0.6 mm
Clear and precise imaging

5 kW
Ready for demanding surgical applications


200 kHU
Anode heat capacity, build for long-lasting performance



Small Pixel Size · High Matrix

Micron Level Pixel Pitch
HD dynamic FPD for clearer image details

9×9/12×12 inch
Detector size for wider FOV and smoother dynamic imaging



8 MP Medical Monitor

I-Dose Low Dose Control

High-DQE Flat Panel Detector

Dedicated Pediatric Low-Dose Mode

Dual-Laser Positioning System

Four-Level Adjustable Dose

Detachable Grid

Anatomy-Specific EPX Curves

Combined Spectral Filtration

SMART Image Processing Algorithm

Adjustable Digital Pulse Rate/Width

Compartment & Wireless Exposure Foot Switch

Integrated X-ray Generator

Real-Time DAP Display

Exposure Time Reminder

