

Veterinary High-power Portable X-ray Generator

TUTOM50E-8A

- **8kW High-frequency generator**
Ensuring reliable and powerful performance every time.
- **7-inch large touch screen**
Streamlining clinical practice and optimizing animal care.
- **Built-in collimator design**
Improved radiation shielding for safety at the interface.
- **Applicable in different scenarios**
Perfect for use in veterinary indoor and outdoor settings.



Veterinary DR Solution

SONTU

Technical Specifications

Items	Parameters
Power Source	Single phase, 220V/110V AC
Input Power	1kVA
Maximum Power	8kW (80mA, 100kV)
Instaneous Power	3kVA
Inverter Frequency	100kHz
Tube Voltage	40kV ~ 125kV
Tube Current	3.2mA - 160mA
Loading Time	1~10000ms
Tube Current-time Product	0.1~250mAs
Focal Spot	2.0mm
Target Anode Angle	15°
Anode Heat Capacity	30kJ (43kHU)



Advantages

Prompt adjustment of exposure parameters and APR setting.

A wireless hand switch that enables remote exposure from up to 15m.



Wireless Hand Switch



Shenzhen SONTU Medical Imaging Equipment Co., Ltd.

Add: Building 1, Shenzhen Biomedicine Innovations Industrial Park, No. 14, Jinhui Road,
Pingshan New District, Shenzhen, 518122, China

Tel: +86 755 27905585

Fax: +86 755 27909350

E-mail: global@sontu.com



SONTU

Medical Imaging